



Magnetic stirrer RCT standard 2.0

Product description	<p>The next generation of our top-selling IKA/LLG-Labware co-branded LLG-RCT Standard magnetic stirrer. New design, new features! Trusted and reliable quality!</p> <ul style="list-style-type: none"> - Digital controlled speed and temperature - Clearly arranged LED display - Digital display of set safety temperature limit - Hot top indicator, hot surface warning to prevent burns - Adjustable safety circuit of heating plate temperature - Digital error code display - Connection for external PT1000 temperature probe. The control accuracy with temperature sensor is ± 1 °C. - Aluminium hotplate for optimal heat transfer - Enhanced safety due to improved heat control technology - USB interface - Powerful motor - Improved magnetic coupling - Supplied with housing protection cover - Tested acc. to DIN EN IEC 61010-1 																						
Scope of supply	Magnetic stirrer, temperature probe PT1000, screwdriver, protective cover, C16 mains cable (H11)																						
Specifications	<table> <tr> <td>Stirring quantity (H₂O):</td><td>20 l</td></tr> <tr> <td>Speed range:</td><td>50 ... 1500 rpm</td></tr> <tr> <td>Heating plate:</td><td>Ø 135 mm</td></tr> <tr> <td>Temperature range:</td><td>RT + device self-heating ... 310 °C</td></tr> <tr> <td>Heat output:</td><td>600 W</td></tr> <tr> <td>Safety temperature range:</td><td>50 ... 360 °C</td></tr> <tr> <td>Dimensions (W x D x H):</td><td>160 x 270 x 85 mm</td></tr> <tr> <td>Weight:</td><td>2.4 kg</td></tr> <tr> <td>Power supply:</td><td>220 ... 230/115/100 V, 50/60 Hz</td></tr> <tr> <td>IP code:</td><td>IP 42</td></tr> <tr> <td>Warranty:</td><td>3 years</td></tr> </table>	Stirring quantity (H₂O):	20 l	Speed range:	50 ... 1500 rpm	Heating plate:	Ø 135 mm	Temperature range:	RT + device self-heating ... 310 °C	Heat output:	600 W	Safety temperature range:	50 ... 360 °C	Dimensions (W x D x H):	160 x 270 x 85 mm	Weight:	2.4 kg	Power supply:	220 ... 230/115/100 V, 50/60 Hz	IP code:	IP 42	Warranty:	3 years
Stirring quantity (H₂O):	20 l																						
Speed range:	50 ... 1500 rpm																						
Heating plate:	Ø 135 mm																						
Temperature range:	RT + device self-heating ... 310 °C																						
Heat output:	600 W																						
Safety temperature range:	50 ... 360 °C																						
Dimensions (W x D x H):	160 x 270 x 85 mm																						
Weight:	2.4 kg																						
Power supply:	220 ... 230/115/100 V, 50/60 Hz																						
IP code:	IP 42																						
Warranty:	3 years																						

Type	PK	Cat. No.
LLG-RCT standard 2.0, EU plug	1	6.320 400
LLG-RCT standard 2.0, CH plug	1	6.320 401
LLG-RCT standard 2.0, UK plug	1	6.320 402