



Pipettierhelfer easy

Produktbeschreibung

Der LLG-Pipettierhelfer easy ermöglicht schnelles und exaktes Pipettieren mit Glas- und Kunststoffpipetten.

- Ergonomische und leichte Konstruktion für konstanten Bedienkomfort bei längerem Gebrauch
- Einzigartiger Winkel des Nasenkonus für mehr Flexibilität im Betrieb, z. B. in Abzugshauben
- Umweltfreundliche, wiederaufladbare NiMH-Batterien, die bis zu 8 Stunden Dauerbetrieb ermöglichen
- Leicht austauschbare Batterien
- Autoklavierbarer Spitzenkonus zur Vermeidung von Kreuzkontaminationen

Lieferumfang

Pipettierhelfer, Handbuch, Ladegerät, Wandhalterung, Tischständer

Typ

LLG-Pipettierhelfer easy, EU-Stecker

LLG-Pipettierhelfer easy, UK-Stecker

VE Bestell Nr.

1 6.263 683

1 6.263 684



Pipette controller easy

Product description

The LLG-Pipette controller easy enables fast and exact pipetting with glass and plastic pipettes.

- Ergonomic and lightweight construction for constant ease of use over prolonged usage
- Unique nose cone angle for enhanced flexibility in operation e.g. in fume hoods
- Eco-friendly rechargeable NiMH batteries, offering up to 8 hours of continuous working
- Easily replaceable batteries
- Autoclavable tip cone for avoiding cross-contamination

Scope of supply

Pipette controller, manual, charger, wall-mount, desktop stand

Type	PK	Cat. No.
LLG-Pipette controller easy, EU plug	1	6.263 683
LLG-Pipette controller easy, UK plug	1	6.263 684