



uniPERISTALTICPUMP 1

Product description	The extra compact peristaltic pump. Ideal for continuous dosing applications. Space-saving, robust and with a wide range of easily exchangeable pump heads.	
Characteristics	<ul style="list-style-type: none"> - Stable flow rates - Compact design - Reversible working direction - Pump heads easily exchangeable - Various pump heads with different hose diameters available as accessory - Low power consumption - Corrosion-resistant ABS housing - OLED display for indication of speed and working direction 	
Delivery includes	Pump, pump head, Pharmed® hose connectors 3 mm ID, 5 mm AD	
Specifications	Speed range:	0.1 ... 150 min ⁻¹
	Speed accuracy:	±0.1 min ⁻¹
	Flow rate range:	0.004 ... 63.69 ml/min
	Weight:	0.8 kg
	Dimensions (W x D x H):	180 x 115 x 125 mm
	IP code:	IP 31
	Working temperature range:	0 ... 40 °C
	Max. humidity:	80 %
	Power rating:	10 W
	Power supply:	100 ... 240 V, 50/60 Hz
	Warranty:	3 years

Type	PK	Cat. No.
LLG-uniPERISTALTICPUMP 1, with EU plug	1	6.263 540
LLG-uniPERISTALTICPUMP 1, with UK plug	1	6.263 541