



uniPERISTALTICPUMP 3

Product description	The peristaltic pump with high pumping capacity. Ideal for continuous dosing applications. Space-saving, robust and with stable flow rate.	
Characteristics	<ul style="list-style-type: none"> - High pumping capacity of 398 ml/min - Reversible working direction - Hose connectors easily exchangeable - Low power consumption - Corrosion-resistant ABS housing - OLED display for indication of speed and working direction 	
Delivery includes	Pump, pump head, Pharmed® hose connectors 6,4 mm ID, 9,6 mm AD	
Specifications	Speed range:	0.1 ... 150 min ⁻¹
	Speed accuracy:	±0.1 min ⁻¹
	Flow rate range:	0.08 ... 398 ml/min
	Weight:	0.8 kg
	Dimensions (W x D x H):	190 x 128 x 135 mm
	IP code:	IP 31
	Working temperature range:	0 ... 40 °C
	Max. humidity:	80 %
	Power rating:	20 W
	Power supply:	100 ... 240 V, 50/60 Hz
	Warranty:	3 years

Type	PK	Cat. No.
LLG-uniPERISTALTICPUMP 3, with EU plug	1	6.263 550
LLG-uniPERISTALTICPUMP 3, with UK plug	1	6.263 551