



uniTEXER 4

Product description	Compact vortexer with brushless motor and infinitely variable speed control up to 4200 rpm for optimal mixing results. The digital display alternately shows the set speed and time. <ul style="list-style-type: none">- Programmable Pulse mode- 3-way switch (on/off/touch)- Timer adjustable from 1 min up to 999 min	
Scope of supply	Vortexer, conical tube attachment, round disk attachment, universal plug	
Specifications	Type of movement:	Orbital
	Orbital diameter:	4 mm
	Speed range:	300 ... 4200 rpm
	Dimensions (W x D x H):	205 x 136 x 138.5 mm
	Weight:	2.8 kg
	Power supply:	110 ... 240 V, 50/60 Hz
	IP code:	IP 42
	Warranty:	3 years

Type

LLG-uniTEXER 4

PK

1

Cat. No.

6.263 500