



LLG-pH Meter 7

Product description	The LLG pH Meter 7 is an easy-to-use, robust table-top instrument for measuring pH and redox potential. The large color-touch-display and the intuitive, language-independent user interface ensure fast and easy use. Particularly suitable for routine applications in laboratories, production plants and schools.	
Characteristics	<ul style="list-style-type: none"> - Available as set with gel combination electrode („Made in Germany“) with temperature sensor - Large color-touch-display with language-independent graphical user interface - Automatic temperature compensation - Memory for 500 values - 3-point calibration - Reminder function for upcoming calibrations - BNC connection for the electrode - RS232 interface for data transfer - Electrode holder can be mounted left and right 	
Delivery includes	LLG pH meter 7: pH meter, electrode holder, protective cover, power supply unit LLG-pH Meter 7 set: pH meter, combination electrode, electrode holder, protective cover, power supply unit	
Specifications	pH range: pH resolution: pH accuracy: Temperature range electrode: Temperature range pH meter: mV range: mV resolution: mV accuracy: Dimensions (W x D x H): Weight: Protection class: Power supply: Warranty:	0 ... 14 0.01 ±0.01 0 ... 80 °C -5 ... 105 °C -2000 ... 2000 mV 1 mV ±1 mV 180 x 240 x 80 mm 0.7 kg IP 54 9 V DC/1A 3 years

Type	PK	Cat. No.
LLG-pH Meter 7, with EU plug	1	6.263 600
LLG-pH Meter 7 set, with EU plug	1	6.263 601
LLG-pH Meter 7, with UK plug	1	6.263 606
LLG-pH Meter 7 set, with UK plug	1	6.263 607
UK plug for LLG-pH Meter 7	1	6.263 634