



uniSTIRRER 2

Product description	Compact magnetic stirrer without heating for volumes up to 1000ml. Adjustable speed from 0 to 2000 rpm. Plate diameter 120 mm. Body made of ABS. Support rod, aluminium (Cat. No. 6.263 435) available as an optional accessory.																								
Specifications	<table> <tr> <td>Stirring quantity (H₂O):</td><td>1000 ml</td></tr> <tr> <td>Speed range:</td><td>0 to 2000 rpm</td></tr> <tr> <td>Permissible ambient temperature:</td><td>0 to 50 °C</td></tr> <tr> <td>Magnetic bar:</td><td>25 x 8 mm</td></tr> <tr> <td>Display:</td><td>Analog</td></tr> <tr> <td>Heating plate:</td><td>120 x 120 mm</td></tr> <tr> <td>Material:</td><td>ABS case</td></tr> <tr> <td>Dimensions (W x D x H):</td><td>150 x 135 x 53 mm</td></tr> <tr> <td>Weight:</td><td>0.5 kg</td></tr> <tr> <td>Power supply:</td><td>100 - 230 V, 50/60 Hz</td></tr> <tr> <td>Protection class:</td><td>IP 21</td></tr> <tr> <td>Warranty:</td><td>3 years</td></tr> </table>	Stirring quantity (H₂O):	1000 ml	Speed range:	0 to 2000 rpm	Permissible ambient temperature:	0 to 50 °C	Magnetic bar:	25 x 8 mm	Display:	Analog	Heating plate:	120 x 120 mm	Material:	ABS case	Dimensions (W x D x H):	150 x 135 x 53 mm	Weight:	0.5 kg	Power supply:	100 - 230 V, 50/60 Hz	Protection class:	IP 21	Warranty:	3 years
Stirring quantity (H₂O):	1000 ml																								
Speed range:	0 to 2000 rpm																								
Permissible ambient temperature:	0 to 50 °C																								
Magnetic bar:	25 x 8 mm																								
Display:	Analog																								
Heating plate:	120 x 120 mm																								
Material:	ABS case																								
Dimensions (W x D x H):	150 x 135 x 53 mm																								
Weight:	0.5 kg																								
Power supply:	100 - 230 V, 50/60 Hz																								
Protection class:	IP 21																								
Warranty:	3 years																								

Description	PK	Cat. No.
LLG-uniSTIRRER 2 with EU-plug	1	6.263 430
LLG-uniSTIRRER 2 with UK-plug	1	6.263 431
Support rod, 25 cm, aluminum	1	6.263 435