



uniSTIRRER 7 pro

Product description	Universal magnetic stirrer with an outstanding price-performance ratio and new design. High level of safety thanks to „Hot“ display with continuing stirring function until the plate temperature has dropped below 50 °C, even when the device is switched off. Integrated control function for precise control of the media temperature by a connected PT1000 temperature probe.	
Characteristics	<ul style="list-style-type: none">-Maintenance-free, brushless DC motor-Stainless steel plate with glass ceramic coating-Heatable-Sturdy metal housing-Powerful stirring-Digital control, backlit display-RS232 for stirrer control and display of all parameters via PC-Display for hot surfaces, even when the device is switched off-Fixed safety circuit of the hotplate temperature at 580 °C-Connection for PT1000 temperature probe, integrated control function-Rod can be mounted on the right and left side of the magnetic stirrer	
Delivery includes	Magnetic stirrer, PT1000 temperature probe, rod, holder, user manual	
Specifications	Stirring quantity (H ₂ O):	20 l
	Speed range:	100 to 1500 rpm
	Accuracy:	±1 rpm
	Magnetic bar:	80 mm
	Display:	LCD
	Heating plate:	184 x 184 mm / stainless steel plate with glass ceramic
	Permissible ambient temperature:	5 to 40 °C
	Temperature range:	Ambient to 550 °C
	Setting accuracy:	±0.1 °C
	Control accuracy:	±0.2 °C
	Heat output:	1000 W
	Safety temperature:	580 °C
	Dimensions (W x D x H):	215 x 360 x 112 mm
	Weight:	5.3 kg
Power supply:	200 - 240 V, 50/60 Hz	
Protection class:	IP 42	
Warranty:	3 years	

Type	PK	Cat. No.
LLG-uniSTIRRER 7 pro, EU plug	1	6.263 448
LLG-uniSTIRRER 7 pro, UK plug	1	6.263 449