



uniSTIRRER 1 and 1/M4

Product description	An ultra-flat inductive magnetic stirrer without heating with timer. The modern magnet coil technology consists only of non-moving parts and therefore guarantees a wear-free drive. For a better mixing of samples it is possible to change the direction of rotation every 30 seconds.	
Characteristics	<ul style="list-style-type: none"> - Timer (1min to 99min, infinitely) - 100% wear- and maintenance-free inductive drive - LED-display for an exact adjustment of the rotation speed - Digital regulated rotation speed - Alteration of the stirring direction - High protection class IP 65 - Stainless steel surface 	
Specifications		
	LLG-uniSTIRRER 1	LLG-uniSTIRRER 1/M4
Stirring positions:	1	4
Stirring quantity (H₂O):	800 ml	4 x 400 ml = 1600 ml
Distance between Stirring positions:	-	95 mm
Speed range:	15 to 1500 rpm	15 to 1200 rpm
Magnetic bar:	max. 25 mm	max. 25 mm
Display:	LED	LED
Heating plate:	94 mm diam./ Stainless steel	(W x D) 122 x 385 mm
Permissible ambient temperature:	5 to 40 °C	5 to 40 °C
Dimensions (W x D x H):	110 x 155 x 12 mm	126 x 457 x 13,5 mm
Weight:	0.5 kg	1.9 kg
Power supply:	100 - 230 V, 50/60 Hz	100 - 240 V, 50/60 Hz
Protection class:	IP 65	IP 65
Warranty:	3 years	3 years

Type	PK	Cat. No.
LLG-uniSTIRRER 1	1	6.263 420
LLG-uniSTIRRER 1/M4, with EU plug	1	6.263 425
LLG-uniSTIRRER 1/M4, with UK plug	1	6.263 426