



uniPLATESHAKER 1

Product description	The LLG-uniPLATESHAKER 1 has been designed for accommodating and mixing a single micro plate. Due to its small footprint it can be used in working areas with limited space such as fume hoods and bio safety cabinets. It can also be used inside a cold-room.																								
Characteristics	<ul style="list-style-type: none"> - Digital Display with Speed and Time settings - Variable speed from 300 to 1800 rpm - Maintenance free brushless DC motor - Orbital Diameter of 2 mm for thorough mixing - Timer setting from 1 to 99 min and continuous mode - Universal fit, comes with interchangeable clamps for different micro plate sizes - Comes with foam attachment for 20 x 1.5/2 ml microtubes 																								
Specifications	<table> <tr> <td>Type of movement:</td><td>orbital</td></tr> <tr> <td>Orbital diameter:</td><td>2 mm</td></tr> <tr> <td>Speed range:</td><td>300 to 1800 rpm</td></tr> <tr> <td>Max. load:</td><td>1 kg</td></tr> <tr> <td>Motor:</td><td>Brushless DC motor</td></tr> <tr> <td>Timer:</td><td>1 to 99 min, continuous</td></tr> <tr> <td>Power consumption:</td><td>15 W</td></tr> <tr> <td>Permissible ambient temperature:</td><td>4 to 40 °C</td></tr> <tr> <td>Dimensions (W x D x H):</td><td>127 x 85 x 86 mm</td></tr> <tr> <td>Weight:</td><td>1.4 kg</td></tr> <tr> <td>Power supply:</td><td>100-230V, 50/60 Hz</td></tr> <tr> <td>Warranty:</td><td>3 years</td></tr> </table>	Type of movement:	orbital	Orbital diameter:	2 mm	Speed range:	300 to 1800 rpm	Max. load:	1 kg	Motor:	Brushless DC motor	Timer:	1 to 99 min, continuous	Power consumption:	15 W	Permissible ambient temperature:	4 to 40 °C	Dimensions (W x D x H):	127 x 85 x 86 mm	Weight:	1.4 kg	Power supply:	100-230V, 50/60 Hz	Warranty:	3 years
Type of movement:	orbital																								
Orbital diameter:	2 mm																								
Speed range:	300 to 1800 rpm																								
Max. load:	1 kg																								
Motor:	Brushless DC motor																								
Timer:	1 to 99 min, continuous																								
Power consumption:	15 W																								
Permissible ambient temperature:	4 to 40 °C																								
Dimensions (W x D x H):	127 x 85 x 86 mm																								
Weight:	1.4 kg																								
Power supply:	100-230V, 50/60 Hz																								
Warranty:	3 years																								

Description

Description	PK	Cat. No.
LLG-uniPLATESHAKER 1 with EU-plug	1	6.263 488
LLG-uniPLATESHAKER 1 with UK-plug	1	6.263 489