



uniTEXER 1

Product description	Compact test tube shaker with brushless motor drive and touch function for short-term operation. Adjustable speed control up to 4500 rpm for excellent mixing results in a very short time. Especially for small volumes.	
Specifications	Type of movement:	orbital
	Orbit:	4.5 mm
	Speed range:	0 bis 4500 rpm
	Permissible ambient temperature:	5 to 40 °C
	Permissible relative moisture:	< 80 %
	Dimensions (W x D x H):	100 x 100 x 70 mm
	Weight:	0.55 kg
	Power supply:	100 - 240 V, 50/60 Hz
	Protection class:	IP43
	Warranty:	3 years

Type	PK	Cat. No.
LLG-uniTEXER 1, with EU-plug	1	6.263 490
LLG-uniTEXER 1, with UK-plug	1	6.263 491
Silicone cap for LLG-uniTEXER 1	1	6.263 492